

WinCare Release Notes

WinCare Version 6.0 Service Release 3 - November 19, 2009

General:

- New MDS and Swing Bed Case Mix Index C05 and C06 that are effective October 2009. "Due to the new CMI change, facilities need to do the following for each Medicare assessment with reference date (A3a) of 10/1/2009 or later and that has RMX Medicare RUGS score:
 1. If this is a comprehensive assessment (assessment with section V), note the VB2
 2. Open the assessment in Section R.
 3. Click the calculate button for T3a. IF this assessment meets the special conditions in the CMI change, then it will change to RVL otherwise will remain unchanged.
 4. If this is a comprehensive assessment, fill in the VB2 with the date noted in Step #1.
 5. Check & Close the assessment in Section R. Note that the Closed date is NOT transmitted to the state, the state considers assessments locked the day they are successfully accepted into the State MDS system.
 6. If the assessment DID change to RVL RUGS score, be sure to inform the Medicare Biller of the corrected RUGS score for this assessment so the biller can put the correct RUGs on the claim.
 7. Repeats steps #1-7 for all assessments that are Medicare, RMX RUGs score, and ARD is 10/1/2009 or later."
- New Texas CMI set for Medicaid.
- Includes the 2010 ICD-9-CM codes that are effective October 2009.
- Updates to the MDS to fix a problem with Section K on certain Other State assessments, and a problem opening section M1 on MPAF assessments.
- An option was added to Notes that allows printing of the Progress Notes by date range.
- A Security Report has been added to the Security Options Screen, that reports the User's Names, Logon Name, Security Level, and Title. To access the report go to Utilities, Facility, WinCare Security, Security Report. Note there is no preview option for the report. If you wish to preview the report you will need a third party application like FinePrint.
- Enhancement to Administration Times: When an Administration Time is changed in the MAR Setup, of Utilities, the change will no longer affect discontinued or invalidated MAR records that are assigned that Administration Time. The change will affect all MAR records with the Administration Time Frequency Code that have not previously been customized. If you wish to change an Administration Time option without affecting existing MAR records, then you must add a new Administration Time Frequency Code.
- The location that the WinCare EXE Information report is saved to, has been changed to %AppData%\Wincare
- Enhancement to the CQI module, sometimes old matrixes didn't show the resident's name when they were discharged. Now discharged residents names always show up.